



MedicalDirector®

Importing a Referral Template into Medical Director

1. Download and save the referral template to an accessible location, such as your Desktop.

IMPORTANT: Do not open the referral template before importing it, as this may corrupt the file.

2. In Medical Director, go to Tools > Letter Writer.
3. Select File > New.
4. Select Blank Template and click OK. A new blank document will open.
5. Select File > Import.
6. Locate and select the file, click Open.
7. Select File > Save as Template. The Save window appears.
8. Enter a name for the template and click Save. The template is now ready for use