

Collodion Adhesive

The adhesive of choice for long-term EEG/EP studies

Collodion is the adhesive of choice for long-term or any EEG/EP studies where secure and stable electrode application is required.

- ✓ Specially formulated by SLE for EEG use
- ✓ Supplied in convenient 6.5 ml individual tubes
- ✓ Cleaner and safer than dispensing from larger containers
- ✓ Supplied in airtight plastic containers
- ✓ Made on site at the SLE factory

Collodion Adhesive is a specially formulated mixture of diethyl ether and methylated spirit for application of EEG electrodes to the head.

It is supplied in 6.5 ml tubes, and is available in airtight plastic tubs of 20 or 50 tubes. If kept in a cool environment Collodion has a shelf life of at least 6 months. SLE Collodion Adhesive is produced at the SLE factory in the United Kingdom.

Collodion LV

A lower viscosity version of Collodion, called Collodion LV is also available. It has a thinner consistency and is available in 6.5 ml tubes in tubs of 20.

Collodion and Collodion LV can be removed with Mavidon Adhesive Remover.

Ordering Details

Part Number	Description
B0005	Collodion adhesive 6.5 ml tube. (Tub of 20)
B0005/50	Collodion adhesive 6.5 ml tube. (Tub of 50)
B0005/09	Mavidon medical adhesive remover for Collodion (473ml)
B0112	Electrode applicator air gun
B0113	Electrode applicator gun tips - each
B0005/0LV/020	Collodion LV adhesive 6.5 ml tube. (Tub of 20)
Stick on Electrodes with 2mm Plug	
B0165/10	Chlorided silver 6mm electrode, 1m cable, 2mm plug - set of 24
B0195/100/024	Chlorided silver 9mm electrode, 1m cable, 2mm plug - set of 24
B0162/50	Unchlorided silver 6mm electrode, 1m cable, 2mm plug - set of 24
B0192/02	Unchlorided silver 9mm electrode, 1m cable, 2mm plug - set of 24
B0192/15	Gold plated silver 9mm electrode, 1m cable, 2mm plug - set of 24
Stick on Disc Electrodes with Touch Proof Plug	
B0166/100/024	Chlorided silver 6mm disc electrode, 1m cable - set of 24
B0196/100/024	Chlorided silver 9mm disc electrode, 1m cable - set of 24
B0166/50/100/024	Unchlorided silver 6mm disc electrode, 1m cable - set of 24
B0196/50/100/024	Unchlorided silver 9mm disc electrode, 1m cable - set of 24
B0167/100/024	Gold plated Silver 6mm disc electrode, 1m cable - set of 24
B0197/100/024	Gold plated silver 9mm disc electrode, 1m cable - set of 24



B0005 - Collodion Adhesive



B0005/0LV/020 - Collodion LV Adhesive



B0005/09 Medical Adhesive Remover

*SLE Collodion is extremely flammable and should be handled with care. Use in a well ventilated area away from sources of ignition.